

Smart Laser & Plasma Systems

32ch- Analog to Digital Converter

2D Time series Temperature and
Concentration distributions measurement system



Model Number

32ch Analog to Digital Converter : AD-32ch-20MHz-12bit-I1M-BNC


Specification

Number of channels	32ch
Vertical resolution	12bits
Maximum sampling rate	20MS/s
Voltage resistance	50Ω
Bandwidth	20MHz
Connector	BNC
Operating voltage	AC 100-240V 50/60Hz
Environment	Temperature : 10~40°C Humidity : 30~80%RH (No Condensation)


Dimensions

● Analog to Digital Converter

Front view



Rear view



Smart Laser & Plasma Systems Co.
VBL, 2nd Floor, 2-1, Minamijyosanjima,
Tokushima, 770-8506 JAPAN

TEL: +81-88-656-7556 **FAX: +81-88-656-7556**
E-mail: info@slps.co.jp
WEB: <http://slps.co.jp/index.html>